

**X-ray Irradiator Comparison Chart - Blood**  
5/15/2018

Carolyn MacKenzie, [cjmackenzie@berkeley.edu](mailto:cjmackenzie@berkeley.edu)  
University of California

Make	Primary Use	Model	Energy	Typical Dose Rates Dose Rates	Key Features	Price Quote Base System	Shipping/ Training and Installation	Other Key Options	Initial Warranty Covered	Warranty 13-24 Mo.	Warranty 24-36 Mo.	Warranty 36-48 Mo.	External Radiation Field at max kVp and mA	Weight
RadSource <sup>1</sup>	Blood	RS 3400	155 Kv / 11 mA / QUASTAR	holds up to 6 (1 liter) blood canisters that hold 1 300 ml blood bag each with a 25 Gy dose in <5 minutes	stand alone-self-contained cooling- plug an play, on wheels, - FDA approved. 1 x-ray tube QUASTAR in center and canisters rotate around it.		shipping	barcode scanner \$4K; syringe holders \$55., extra canisters \$1.5k, phantoms	1 yr				<0.5 mR/hr at 5 sm off external surfaces and <0.1 mR/hr at 10 cm	2450 lbs
Best Theratronics <sup>2</sup>	Blood	MK1	160 kV	delivers 25 Gy in 6 minutes to a 1.7L canister so it irradiates 3 red cell bags in 4 minutes. 6.3 Gy/min	stand alone- 1 x-ray tube; self-contained cooling system		Install/Training included / Shipping included	optional bar code reader \$6.3K;	1 yr				<0.5 mR/hr at 5 sm off external surfaces and <0.1 mR/hr at 10 cm	1309 lbs
Best Theratronics <sup>2</sup>	Blood	MK2	160 kV	lg canisters 2.0 l (holds 4 300 ml bags) or 3.5 l size (holds 6-8 300 ml bags) - delivers 25 Gy in 3.8-5.7 minutes depending on canister size; 5-9 Gy minute depending on the canister size	stand alone- has 2 x-ray tubes above and below the sample for uniform dose.		Install/Training included / Shipping included	database management system \$10.6K; heat exchanger for had water \$16.1K; optional water chiller \$22.3K	1 yr				<0.5 mR/hr at 5 sm off external surfaces and <0.1 mR/hr at 10 cm	2200 lbs

**1. RAD SOURCE TECHNOLOGIES** 480 BROGDON RD. SUITE 500 SUWANEE, GA 30024  
E-MAIL: [RADAMS@RADSOURCE.COM](mailto:RADAMS@RADSOURCE.COM) Mr. Richard Adams  
PHONE: 815-477-1291  
WEBSITE: [HTTP://WWW.RADSOURCE.COM/](http://WWW.RADSOURCE.COM/)

**2. BEST THERATRONICS** 413 MARCH RD OTTOWA, ONTARIO K2K 0E4 CANDADA  
E-MAIL: [DAN.AITKENHEAD@THERATRONICS.CA](mailto:DAN.AITKENHEAD@THERATRONICS.CA) (Mr. Dan Aitkenhead) PHONE: 613-591-2100  
WEBSITE: [HTTP://WWW.THERATRONICS.CA/](http://WWW.THERATRONICS.CA/)